


# Potentiometric Electrodes: Glass pH Electrode

Gabriel Ortega, PhD

[gortega@chem.ucsb.edu](mailto:gortega@chem.ucsb.edu)


University of Toronto Mississauga

Dec 4<sup>th</sup> 2020



Why is pH important?

## Biology



Blood: pH  $\sim$  7.35 – 7.45  
Tumor environment: pH  $<$  7  
Gastric acids: pH  $\sim$  2

## Ecology




Ocean acidification  
Coral bleaching  
Acid rain

## Industry





Food industry  
Fermentation  
Dyes

## Geology



Soil analysis, planetary exploration

## Synthesis



Catalysis and control of reactions

$$\text{pH} = -\log [\text{H}^+]$$

How can we measure pH?

How can we measure the ion concentrations?

# Chemical Potential

rate of change in **energy** of the system with respect to the change in the concentration of species in a chemical reaction

$$\text{chemical potential} = \frac{\text{potential energy}}{\text{number of molecules}}$$

(at constant  $T, P, N$ , partial molar Gibbs free energy,  $\bar{G}$ , in J/mol)

If molecules are charged  $\longrightarrow$  **Electric Potential**

$$\text{electric potential} = \frac{\text{electrostatic potential energy}}{\text{charge}}$$

E (J/C = V)

$$G = G^0 - RT \ln \frac{a_{products}}{a_{reactants}}$$

$$\updownarrow E = G / nF$$

## Nernst Equation

Driving force bringing the reaction to equilibrium

$$E = E^0 - \frac{RT}{nF} \ln \frac{a_{products}}{a_{reactants}}$$

Reaction product similar to an equilibrium constant... for a system that is NOT at equilibrium!

At equilibrium, potential is...

**ZERO**

$E^0$  = standard reduction potential (when  $a_i = 1$ )

R = gas constant,  $8.314 \text{ J K}^{-1} \text{ mol}^{-1}$

T = temperature, 298K

n = moles of electrons exchanged

F = Faraday constant,  $F = q_{e^-} N_A = 96485.33 \text{ C mol}^{-1}$

$a_i$  = activities ( $\approx$  concentrations, M)


How do we measure E?

# Potentiometry

Measurement of a difference in electric potential (a voltage) to obtain chemical information (concentration of charged species)



# Measuring potential differences




$$\Delta E = E_{\text{ind}} - E_{\text{ref}}$$

Ag|AgCl (sat. KCl):  $E_{\text{ref}} = + 0.197 \text{ V}$

$E^0 (\text{Ag}^+ | \text{Ag}^0) = + 0.799 \text{ V}$


$$E = E^0 - \frac{RT}{nF} \ln \frac{a_{\text{products}}}{a_{\text{reactants}}}$$



$$\Delta E = E_{\text{ind}} - E_{\text{ref}}$$


Ag|AgCl (sat. KCl):  $E_{\text{ref}} = + 0.197 \text{ V}$

$E^0 (\text{Ag}^+ | \text{Ag}^0) = + 0.799 \text{ V}$



not generalizable to all metals, unspecific

$$E = E^0 - \frac{RT}{nF} \ln \frac{a_{\text{products}}}{a_{\text{reactants}}}$$



$$E = E^0 - \left[ \frac{8.314 \cdot 298}{n \cdot 96485.33} \cdot 2.303 \right] \log \frac{a_{\text{products}}}{a_{\text{reactants}}}$$

$$E = E^0 - \frac{0.0592}{n} \log \frac{a_{\text{products}}}{a_{\text{reactants}}}$$

$$E = 0.799 - \frac{0.0592}{1} \log \frac{1}{[\text{Ag}^+]}$$


$$E = 0.799 - \frac{0.0592}{1} \log \frac{1}{0.1}$$

$$E_{\text{ind}} = 0.799 - 0.059 \cdot 1 = 0.740 \text{ V}$$

$$\Delta E = 0.740 - 0.197 = 0.543 \text{ V}$$

How Else Can We Generate Potentials?

# Junction Potentials



$\Delta C \rightarrow$  chemical potential

If molecules are charged:

$\Delta C \rightarrow$  electric potential,  $\Delta E$  (V)

*(as opposed to, e.g., fuel cells and batteries, this potential does not involve a redox process, but **mass transport** of charged molecules)*

If the junction is selectively permeable to only one ion





**Ion-Selective Electrode!**

Which Type of Membrane Is Selective for  $H^+$ ?

Glass!!

# Glass pH Electrode




Potential at the membrane

$$E = E^0 - \frac{RT}{nF} 2.3 \log \frac{a_{H^+, sol}}{a_{H^+, ref}}$$

$$E = E^0 - 0.0592 \Delta pH$$

How can we measure this potential?

# Glass pH Electrode



$$\Delta E = E_{ind} - E_{ref}$$

$$\Delta E = E_{ref}^{ind} + E_{membrane} - E_{ref}$$

$$E_{membrane} = 0.0592 (pH_{sol} - pH_{ref})$$

$$E = constant + 0.0592 \cdot pH_{sol}$$

**Beware of junction potentials!**


$$E = constant + 0.0592 \cdot pH_{sol} + E_{junction}$$

**Unknown, not constant**

$$6 \text{ mV} \approx 0.1 \text{ pH units}$$



# Glass pH Electrode



2-3 point calibration in the range of pH of your sample

$$E = E_{ref}^{ind} + E_{membrane} - E_{ref}^{ref}$$

$$E_{membrane} = 0.0592 (pH_{sol} - pH_{ref})$$

$$E = constant + 0.0592 \cdot pH_{sol}$$

**Beware of junction potentials!**

$$E = constant + 0.0592 \cdot pH_{sol} + E_{junction}$$

Unknown, not constant


6 mV  $\approx$  0.1 pH units

# Sources of Error

Calibration solution should be comparable to your sample  
(pH, T, ionic strength)


Selectivity:  $\text{Na}^+$

Acidic error:  
Saturation of glass membrane



Basic error  
 $[\text{H}^+]$  too low  
 $[\text{Na}^+]$  too high

# Glass pH Electrode



# Glass pH Electrode





Phyllis Margaret Tookey Kerridge (1901-1940)

## XCI. THE USE OF THE GLASS ELECTRODE IN BIOCHEMISTRY.

BY PHYLLIS TOOKEY KERRIDGE.

*From the Institute of Physiology, University College, London.*

*(Received June 18th, 1925.)*




Arnold Beckman (1900-2004)

First “portable” pH meter, 1934

Beckman Model G pH Meter



# Glass pH Electrode



# Practical Considerations

Calibration!

Never let dry

To rinse it, gently pat the surface with a wipe.  
Rubbing it will cause static electricity build-up.

The tip is extremely fragile (100  $\mu\text{m}$  thick glass!)

Know when to you need a pH meter:



# References & Links

**Quantitative Chemical Analysis**, 8<sup>th</sup> Ed. Daniel C. Harris

Chapters 13, 14

**Fundamentals of Analytical Chemistry**, 9<sup>th</sup> Ed. Skoog, West, Holler & Crouch

Chapter 21

**Electrochemical Methods: Fundamentals & Applications**, 2<sup>nd</sup> Ed. Bard & Faulkner

Chapter 2 (more advanced, rigorous treatment of electrostatic potentials and junctions)

\*\*\*Bonus Track: Online Tutorial for DIY pH Electrode

<https://www.instructables.com/cheap-DIY-electronic-pH-meter/>



